Hydrogenations with the bpc 2 precise, reproducible, reliable, safe

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Facts and figures

- Dosing of (hydrogen) gas at 1 – 140 bar
- Control mode:
 - constant reaction pressure
 - constant flow rate
- Continuous measuring of H₂ uptake
- Gas consumption in mol, ml
- Accuracy > 99%
- Flowrates 0 – 225 liters n / min
- Functional for:
 - H₂ (Hydrogen)
 - C₂H₄ (Ethylene)
 - CO₂ (Carbon dioxide)
 - CO (Carbon monoxide)
 - C₃H₆ (Propylene)

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The bpc 2 and Büchi pressure reactors – the ultimate solution for hydrogenation

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Built on the experience and the outstanding performance of the first generation hydrogenation system, the bpc 2 offers unmatched reproducibility, proven accuracy, an intuitive interface, and operates at the highest level of safety.

Features

- Reactor (reaction) pressure kept stable during entire experiment
- Extremely precise with accurate reproducibility using our unique volumetric dosing
- Continuous logging of consumed gas (to PC, USB memory stick)
- Touch-screen operation, data display
- Various internal safety interlocks
- Automatic switch-off in event of leakage
- No need for on-site calibration
- Sequence program, automation

Easy operation

1. Purging with inert / active gas
2. Set regulated pressure > delivery pressure
3. Set delivery pressure and start

Simple installation

- 2 lines-in (Gas-supply)
- 1 line-out to reactor
- 1 vent line

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
Setup

- bpc 2 controller**
- Touch-screen operation
 - USB data memory port
 - RS232 port (for optional bls software)
 - Ethernet

- bpc 2 dosing unit**
- Volumetric gas dosing unit (VD*)
 - Pressure control devices
 - Safety interlock valves

Options

- Measure / display of temperature, pressure, stirrer speed, alarm functions
- Sequence program (venting valve included): Example for system preparation before run:
 1. N2 purging
 2. Leak Test
 3. H2 purging
- Bls software for bpc 2 remote control and data storage from and to PC

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
Easy and safe operation – valuable information

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
The Büchi hydrogenation solution

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The bpc 2 is suitable for Büchi pressure reactors of different volumes, pressure ranges and materials.

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Applications

- Hydrogenations (catalytic)
- Chemical research
- Process development
- Scale-up

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* delivery pressure is equal to reaction pressure / *delivery pressure entspricht Reaktionsdruck*

Hydrogenation in glass pressure reactor /
Hydrierung in Glasdruckreaktor

Büchi AG
Gschwaderstrasse 12
CH-8610 Uster/ Switzerland

Phone +41 (0) 44 905 51 11
Fax +41 (0) 44 905 51 22
www.buchiglas.com

 **büchiglasuster**[®]
switzerland Pilot Plant and Reactor Systems